



EASE OF USE, FLEXIBILITY, ERGONOMICS, ROBUSTNESS AND HIGH DENSITY IN A COMPACT FOOTPRINT: THESE ARE THE KEY FACTORS THAT DRIVE OUR DESIGN TO PROVIDE CUSTOMERS WITH THE HANDIEST SOLUTIONS

IWT IS THE IDEAL PARTNER FOR THE CARE OF YOUR FACILITY AND FOR THE OPTIMIZATION OF YOUR DAILY WORKLOAD.

FOR MORE THAN 20 YEARS IWT HAS BEEN DESIGNING, DEVELOPING AND ENGINEERING SOLUTIONS FOR THE LABORATORY ANIMAL INDUSTRY. TODAY, WITH MILLIONS OF ITEM HANDLED EVERY DAY IN THE MOST PRESTIGIOUS VIVARIUMS WORLDWIDE, IWT IS PROUD TO PROVIDE YOU WITH THE MOST COMPREHENSIVE PORTFOLIO OF LOGISTICS SOLUTIONS IN THE MARKET. OUR TEAM OF SPECIALISTS AND ENGINEERS HAS DEVELOPED CRATES, TROLLEYS, RACKS AND ACCESSORIES WITH OPERATORS' CARE AS CENTRAL FULCRUM.

TRANSPORTING AND STOCKING, WASHING AND AUTOCLAVING ARE JUST FEW OF THE TYPICAL STAGES IN THE CLOSE LOOP OF ANY ANIMAL FACILITY. THE ACCURATE AND PROPER SELECTION OF THE RIGHT LOGISTICS IS CRUCIAL TO GUARANTEE A SMOOTH TURNAROUND OF GOODS AVOIDING UNPLEASANT DOUBLE HANDLINGS.

YOUR LOCAL REPRESENTATIVE WILL LEAD YOU TOWARDS THE CHOICE OF THE BEST STANDARD SOLUTION OR TAILOR-MADE CUSTOMIZATION TO FIT YOUR SPECIFIC NEEDS.

INDEX



5
BOTTLE
HANDLING

CRATES
TROLLEYS
ACCESSORIES



9
CAGE
HANDLING

TROLLEYS
RACKS



13
PRESENTATION
RACKS

CAGES
TOPS & WIRELIDS
BOTTLES




16

PRODUCTS
PROFILE



BOTTLE HANDLING



CRATES



TROLLEYS



ACCESSORIES

BOTTLES REQUIRE CONSTANT
TURNOVER. FRIENDLY LOGISTICS IS A
MUST TO ENSURE A SMOOTH PROCESS.

MODERN FACILITIES NEED EQUIPMENT
CAPABLE TO COVER MULTIPURPOSES
FROM HOLDING ROOMS TO CAGE
WASH AREA ENSURING EASE OF USE
FOR THE OPERATOR AND HIGH DENSITY
FOR STOCKING PURPOSES.





TAKE A LOOK AT IWT WASHING BROCHURE TO FIND OUT OUR BOTTLE WASHING SOLUTIONS




CRATES, TROLLEYS AND ACCESSORIES

- Drawer-style transport trolley (8, 10 and 12 spots) for 18 positions bottle crates
- Easy to manoeuvre thanks to their light design
- Crates with embedded sliding lid featuring quick lock/unlock system for easy access to individual bottles
- Different crate sizes to hold from 260 to 700ml (from 9.8 to 23.7oz)
- A set of miscellaneous crates is available for caps and reusable enrichments handling

[9LBCBM18263](#)



[9LBCA003](#)




MULTIPLE SIZES MISCELLANEOUS CRATE



Drawer-style
trolleys
+
Bottle crates
with embedded
sliding lid
=
no more
heavy liftings!

- Trolleys and crates described above:
 - avoid any heavy lifting during change out operations in the holding rooms
 - fit the full facility loop being compatible with manual and automated washing systems
- Washable and autoclavable covers (10µm filtration efficiency):
 - load protection during stocking
 - front access to avoid cover removal during usage



ACCACBO1036

EMBEDDED LOCK/UNLOCK BOTTLE HOLDING SYSTEM

EASY ACCESS WASHABLE AND AUTOCLAVABLE COVERS



CAGE

HANDLING



TROLLEYS



RACKS


MULTIPURPOSE TROLLEYS OR INDIVIDUAL CAGE TYPE RACKS ARE SPECIFICALLY DESIGNED TO INCREASE STOCK DENSITY AND MAXIMIZE BULK AUTOCLAVES LOAD AVOIDING DOUBLE HANDLINGS.

FLEXIBILITY AND ACCESSIBILITY ARE MEETING EASE OF USE AND LIGHTNESS.

CAGE HANDLING



ACCATMWM02



9CASUTRMT1

ACCATXSUITE

MULTIPLE CAGE HANDLING SOLUTIONS

ACCATUSTT6

ACCATM01



TROLLEYS AND RACKS

- Different designs for different purposes:
 - Multishelves solutions to enhance the flexibility of the load
 - Trolleys for nested items suitable for partial cage change operations
 - Racks for single-type cage suitable for integral cage change operations
- Easy to manoeuvre thanks to their light design
- Compatible with the full loop from cage change operations in the holding rooms through washing and autoclaving
- Ratio between floor space and cage numbers optimized to increase density granting load stability
- Standard design for manual applications and custom solutions for automation or special requirements
- Washable and autoclavable covers (10µm filtration efficiency):
 - allergens containment
 - load protection during stocking
 - front access to avoid cover removal during usage
- Alternative trolleys available for rabbit and guinea pig trays and cages



ACCATG80GM



PRESENTATION RACKS



FOR
CAGES



FOR TOPS
& WIRELIDS



FOR
BOTTLES

LOGISTICS THOUGHT TO REMAIN
IN THE INNER LOOP OF THE CAGE
WASH AREA AND TO WORK IN
CONJUNCTION WITH RACK WASHERS.

RACKS ENGINEERED TO OPTIMALLY
COMBINE PROPER LOAD EXPOSURE
TO WATER GUARANTEEING
EFFICIENT DRIPPING AND
MAXIMIZING THROUGHPUT.



9ACCALDTGM

Can you "fill" it?
IWT robots can!



PRESENTATION RACKS FOR CAGES, TOPS & WIRELIDS AND BOTTLES


- Universal racks to accommodate multi-type of cages enhancing flexibility
- High density design to maximize washing throughput
- Bottle washing kit featuring quick-lock connection and rotating arms
- Available in "standard" or "baby" sizes to suit available operational spaces and meeting operators' needs
- Alternative racks available for rabbit and guinea pig trays and cages
- Presentation rack for aquatic tanks to work in conjunction with the related rack washer option
- Standard design for manual applications and custom solutions for automation or special requirements
- Specific racks for IVC tops washing to increase performance avoiding unpleasant water pooling
- Miscellaneous accessories for special applications to be integrated or attached to the racks above (e.g. AWS and IVC valves washing kits)

TAKE A LOOK AT IWT
AUTOMATION BROCHURE
TO FIND OUT MORE




ACCALS02AA

9CASUPRTLY



BABY SIZE RACKS

ACCALDD



ACCALATLS02HD

STANDARD SIZE RACKS



JWI
TECHNICAL




ACCALBCPR

PRODUCTS PROFILE

BOTTLE CRATES

9LBCBMXXXXX

PURPOSE - 18 Positions bottle crates, suitable for Tspa bottles. Inclusive of embedded lid (lock-unlock sliding patented system)



BC-Matic code	Tecniplast bottle
9LBCBM18263	ACBT0202 - ACBT0252 - ACBT0262 - ACBT0312
9LBCBM18403	ACBT0402
9LBCBM18503	ACBT0502 - ACBT0702
9LBCBM184234	ACBT0312 - ACBT0412

BOTTLE TROLLEYS

ACCACBOXXXB

PURPOSE - Bottle transport trolleys for 18 positions bottle crates. Featuring drawer style crates handling.

OPTIONS - ACCACBOBUM: set of bumpers




1250 x 500 x 1850 mm
49.2 x 19.7 x 72.8 "

CRATE CODE	9LBCBM18263	9LBCBM18403	9LBCBM18503	9LBCBM184234	Cover code
ACCACBO1236	✓	✓	x	✓	ACCACBC02
ACCACBO1036	✓	✓	x	✓	ACCACBC01
ACCACBO0836	x	x	✓	x	ACCACBC01
ACCACBO0836L	✓	✓	x	✓	-


ACCESSORIES

CAPS CRATES

PURPOSE - Sipper tubes crates for washing.



9LBCMBASKS



9LBCBM18CAPS




9LBCMMAGIC

ACCESSORIES


9LBCA00X

PURPOSE - Miscellaneous crates with hinged lid for washing durable environmental enrichments. Suitable also for sipper tubes.



9LBCA003

500 x 500 x 266 mm
19.7 x 19.7 x 10.4 "



9LBCA001

500 x 500 x 123 mm
19.7 x 19.7 x 4.8 "

9CASUTTRMTXX

PURPOSE - Cage bases transport trolley

OPTIONS - 9CASUTTRMT1L: allergen lind for 9CASUTTRMT1

9CASUTTRMT2L: allergen lind for 9CASUTTRMT2

	9CASUTTRMT1		9CASUTTRMT2			
	GM500	GR900	1284L	1285L	1291H	1500U
Cage base	60	36	60	60	45	18
Dimensions	860 x 500 x 280 mm 33.8 x 19.7 x 11"		840 x 500 x 290 mm 33 x 19.7 x 11.4"			

CAGE TROLLEYS



ACCATM01

PURPOSE - Miscellaneous transport trolley.

OPTIONS - ACCATM01R: additional shelf

ACCATBUM: set of bumpers

	1145T	1285L	1291H	1500U	GM500	GR900	ISOCAGE
Cage base	153	186	85	27	138	151	168
Top	232	155	136	75	294	141	120

CAGE TROLLEYS



1250 x 500 x 1858 mm
49.2 x 19.7 x 73.1"

ACCATG02F

PURPOSE - Miscellaneous transport trolley. Mainly designed for integral cage set-ups but flexible to accommodate miscellaneous items.

OPTIONS - ACCATBUM: set of bumpers

	1145T	1285L	1291H	1500U	GM500	GR900	ISOCAGE
Cage base	171	180	64	30	132	51	162
Top	248	155	136	81	294	147	150
Integral cage	64	48	20	18	54	18	48

CAGE TROLLEYS



1250 x 500 x 1858 mm
49.2 x 19.7 x 73.1"

ACCATMWM02

PURPOSE - Miscellaneous transport trolley.

OPTIONS - ACCATMWM02D: doors, suggested in case of integral cages

ACCATMWM02S: intermediate shelf, to split the load for better stability

ACCATMC01: polyester cover suggested in case of open cages transportation

	1145T	1285L	1291H	1500U	GM500	GR900	ISOCAGE
Cage base	126	192	72	20	132	36	180
Top	192	128	108	56	300	80	104
Integral cage	56	48	24	12	54	12	48


CAGE TROLLEYS



1040 x 750 x 1802 mm
40.9 x 29.5 x 70.9"

CAGE TROLLEYS

ACCATUSTTX



1105 x 805 x 1672 mm
43.5 x 31.6 x 65.8"

PURPOSE - Integral cage change procedures. Flexible transport trolley to accommodate different cage types over 6 or 7 shelves.

	ACCATUSTT6			ACCATUSTT7	
	1285L	GM500	GR900	1285L	GM500
Integral cage	48	60	24	56	70

CAGE TROLLEYS

ACCATXSUITE



PURPOSE - Cage bases transport trolley rabbit applications.

OPTIONS - ACCATBUM: set of bumpers

	ACCATXSUITE - X-type	ACCATRSUITE - R-type	ACCATPSUITE - P-type
Cage	6	6	6
Tray	12	12	12
Cage + Tray	4 + 4	4 + 4	4 + 4
Dimensions	738x763x1455mm 29x30x57.2"	765x792x1455 mm 30x31x57.2"	915x807x1455mm 36x31.7x57.2"

CAGE RACKS

ACCATGXXXX



PURPOSE - Integral cage change procedures. Transport trolley dedicated to the specific cage type.

	ACCATG8001		ACCATG80GM	ACCATG42GR	ACCATG16DD
	1284L	1285L	GM500	GR900	GR1800DD
Integral cage	80	80	80	42	16
Dimensions	800 x 1230 x 1900 mm 31.5 x 48.4 x 74.8"		800 x 1105 x 1900 mm 31.5 x 43.5 x 74.8"	800 x 1100 x 1900 mm 31.5 x 43.3 x 74.8"	1006 x 870 x 2033 mm 39.6 x 34.2 x 80"

PRESENTATION RACK FOR TOPS & WIRED LIDS

ACCALPRM 9CASUPRTLY



PURPOSE - "Baby" presentation rack dedicated to IVC tops washing. Two of these are fitting one RW load.

	ACCALPRM				9CASUPRTLY
	1285L	1291H	GM500	GR900	GM500
Integral cage	40	40	40	40	84

PURPOSE - Presentation rack for rabbit and guinea pig applications

ACCALG03	X-type	R-type	P-type	F-SUITE	GP-SUITE
Cage	8	8	8	8	8
Tray	8	8	8	-	8

ACCALG03AA	X-type	R-type	P-type	F-SUITE	GP-SUITE
Cage	4	4	4	4	4
Tray	4	4	4	-	4



ACCALDD

PURPOSE - Presentation rack for GR1800DD applications.

OPTIONS - ACCALDDAFH: crate for food hoppers
 ACCALDDAH: crate for hatches

	GM500	GR900	GR1800DD
Cage base	64	32	32
Top	64	32	32
Sky deck	-	-	32



ACCALXXS02XX

PURPOSE - Universal presentation rack for RW applications.

ACCALS02: 800 series and 900GP

ACCALATLS02: Atlantis

ACCALS02AA: "baby" presentation rack for 800 series, 900GP and Atlantis

ACCALS02	1145T	1285L	1291H	1500U	GM500	GR900	ISOCAGE
Cage base	96	72	56	32	80	40	80
Top	96	72	56	32	80	40	72


ACCALATLS02	1145T	1285L	1291H	1500U	GM500	GR900	ISOCAGE
Cage base	104	80	64	32	88	48	88
Top	104	80	64	32	88	48	88


ACCALS02AA	1145T	1285L	1291H	1500U	GM500	GR900	ISOCAGE
Cage base	48	32	24	8	32	16	32
Top	48	32	24	8	32	16	32



ACCALS02HD	1145T	1285L	GM500	ISOCAGE
Cage base	120	90	90	90
Top	120	90	90	90

ACCALS02HD	1145T	1285L	GM500	ISOCAGE
Cage base	130	100	100	100
Top	130	100	100	100





HEADQUARTERS IWT s.r.l. | Tel. +39 0332.96701 · www.iwtsrl.it · E-mail: info@iwtsrl.it

Sister Companies

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 · www.tecniplast.it · E-mail: info@tecniplast.com.au

CHINA | Tel. + 86 (0)21 50810920 · www.tecniplast.cn · E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. +33 (0)4 72 52 94 41 · www.tecniplast.fr · E-mail: info@tecniplast.fr

GERMANY | Tel. +49 (0)8805 921320 · www.tecniplast.de · E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. + 44 (0)845 0504556 · www.tecniplast.it · E-mail: info@tecniplastuk.com

JAPAN | Tel. +81 (0)3 5770 5375 · www.tecniplast.it · E-mail: info@tecniplastjapan.co.jp

USA | Toll Free: +1 877.669.2243 · www.tecniplastusa.com · E-mail: info@tecniplastusa.com

CANADA | Toll Free: +1 855.347.8718 · www.tecniplastcanada.com · E-mail: info@tecniplastcanada.com

TO FIND YOUR LOCAL DISTRIBUTOR, PLEASE VISIT WWW.IWTSRL.IT

